

Renfert Polish


Diamond polishing pastes



Overview of Renfert Polishing Pastes with recommended polishing tools

Item			Recommended Polishing Tool		
Polishing paste	with dia- mond grit	Material/ Range of Application	Pre-polish	High luster polish	Final Polish
Renfert Polish hybrid materials No. 5103000	✓	Hybrid materials, hybrid ceramic			
Renfert Polish all-in-one No. 5100100	✓	Lithium disilicate, zirconium dioxide, veneering and hybrid ceramics, high performance polymers, veneering composites			
Renfert Polish ZrO ₂ No. 5101100	✓	Zirconium dioxide			
Renfert Polish LiSi ₂ No. 5102100	✓	Lithium disilicate			
Kohinoor L No. 5160001	✓	Veneering ceramic			
Brinell L No. 5190001	✓	Precious PFM alloys, ceramic		*	
Opal L No. 5200001		Veneering composites			
Saphir No. 5150000		Precious metal alloys			
Saphir No. 5150000		Cobalt-chrome alloys			
Universal polishing paste No. 5131000		PMMA, resins			

*Included in Brinell L No. 5190001 scope of delivery



NEW Renfert Polish hybrid materials

Diamond polishing paste

Diamond polishing paste specifically for high-gloss polishing of hybrid materials, for extra-oral use.



Renfert Polish hybrid materials

- Reduced plaque build-up and longer service life of the tooth restorations thanks to optimal surface finish.
- Optimal polishing result thanks to diamond particles in the polishing paste that are specially matched to hybrid materials.

Ordering information

Renfert Polish hybrid materials, 13 g	No. 510 3000
Renfert Polish all-in-one, 18 g	No. 510 0100
Renfert Polish ZrO ₂ , 18 g	No. 510 1100
Renfert Polish LiSi ₂ , 18 g	No. 510 2100

Renfert Polish all-in-one

Universal polishing paste with a wide spectrum of diamond grain sizes, for extra-oral use. Quick, high gloss finish on various materials such as lithium disilicate, zirconium dioxide, veneering and hybrid ceramics such as high performance polymers and veneering composites.

Renfert Polish ZrO₂

Extra-oral polishing paste for zirconium dioxide with specially coordinated abrasive properties due to extra coarse diamond grit.

Renfert Polish LiSi₂

Extra-oral polishing paste specially designed for lithium disilicate with a precise amount of coarse diamond grit.

